beurer

BF 1000 Super Precision

DE	Diagnosewaage Gebrauchsanweisung	2
EN	Diagnostic bathroom scale Instructions for use	20
FR	Pèse-personne impédancemètre Mode d'emploi	37
ES	Báscula diagnóstica Instrucciones de uso	55
IT	Bilancia diagnostica Istruzioni per l'uso	73
TR	Diyanozoz terazisi Kullanım kılavuzu	91
RU	Диагностические весы Инструкция по применению	108
PL	Waga diagnostyczna Instrukcja obsługi	127

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Read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Dear customer,

Thank you for choosing a product from our range. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage, beauty, baby and air.

With kind regards,
Your Beurer team

Getting to know your device

Functions of the device

The purpose of this digital scale is weighing and calculating your personal fitness data.

It is intended for private use.

The scale is equipped with the following functions, which can be used by up to 10 different people:

- Body weight measurement
- Body fat percentage calculation
- Visceral fat percentage calculation (organ fat/abdominal fat)
- Body water content calculation
- Muscle percentage calculation
- Bone mass calculation
- Basal metabolic rate and active metabolic rate calculation.

The scale also includes the following additional features:

- Cross measurement with 8 electrodes allows for segment analysis, i.e. a separate display of individual body parts (arms, legs, torso)
- Innovative menu navigation on large display (in 5 languages),
- Switch between kilograms “kg”, pounds “lb” and stones “st”
- Battery replacement display for low batteries
- Automatic user recognition
- Storage of the last 30 measurements for 10 users if the data cannot be transferred to the app
- Guest mode (without memory function).

The measurement principle

This scale works on the principle of B.I.A. (bioelectrical impedance analysis). This involves the calculation of body content in a matter of seconds, using a current that cannot be felt, is completely harmless and does not pose any risk.

When this measurement of the electrical resistance (impedance) is considered alongside constants and/or individual values (age, height, gender, activity level), the body fat percentage and other variables in the body can be calculated. Muscle tissue and water conduct electricity well and therefore have a lower resistance. In contrast, bones and fat tissue have low conductivity as the fat cells and bones are barely able to conduct the current due to very high levels of resistance.

The electrical contact to the body is enhanced through 4 hand electrodes and 4 foot electrodes to achieve a higher measurement accuracy and better analysis. This sensitive analysis technology means readings may change slightly when measurements are repeated. This is caused by changes during the impedance measurement (e.g. different electrical contact to the electrodes by the hands and feet, different distribution of water in the body).

Please note that the values calculated by the diagnostic bathroom scale only represent an approximation of the actual medical analytical values for the body. Only a specialist doctor can precisely establish the body fat, body water, muscle percentage and bone structure using medical methods (e.g. computer tomography).

Contents

1. Included in delivery.....	21	9. Analysis of measured values	28
2. Signs and symbols	21	10. Evaluating the results	31
3. Warnings and safety notes.....	22	11. Transferring measurements via app	33
4. Device description.....	24	12. Additional information	33
5. Initial use	25	13. Cleaning and maintaining the device	34
6. Configuring the diagnostic bathroom scale with an app.....	26	14. Disposal.....	34
7. Alternatively: configuring the diagnostic bathroom scale without an app.....	26	15. What if there are problems?	34
8. Carrying out the measurement.....	27	16. Technical specifications	35
		17. Warranty / service.....	36



1. Included in delivery

- Diagnostic bathroom scale with raised display on aluminium stand
- 4 x 1.5V type AA batteries
- These instructions for use
- 2 x quick guide
- Mounting kit (4 screws with 2 replacement screws + Allen key)

2. Signs and symbols

The following symbols are used on the device, in these instructions for use, on the packaging and on the type plate for the device:

	WARNING Warning notice indicating a risk of injury or damage to health.
	IMPORTANT Safety note indicating possible damage to the device/accessory.
	Note Note on important information
	Observe the instructions for use
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Dispose of packaging in an environmentally friendly manner
	This product satisfies the requirements of the applicable European and national directives.
	The device may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.
	Manufacturer
	Certification symbol for products that are exported to the Russian Federation and members of the CIS

	<p>Do not load the scale beyond 200 kg/440 lb/31 st.</p>		<p>Danger of slipping: do not step on the scale with wet feet.</p>
	<p>Danger of tipping: position yourself in the centre of the weighing surface.</p>		<p>Place the scale on an even surface. No carpet.</p>

3. Warnings and safety notes

WARNING

- **The scale may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.**
- Do not use during pregnancy.
- Do not step onto the outer edge of the scale on just one side: danger of tipping!
- Swallowing batteries can be extremely dangerous. Keep the batteries and scale out of the reach of small children. Should a battery be swallowed, seek medical assistance immediately.
- Keep packaging material away from children (risk of suffocation).
- Important: do not step onto the scale with wet feet or if the surface of the scale is damp – danger of slipping!

Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
- **⚠ Choking hazard!** Small children may swallow and choke on batteries. Therefore, store batteries out of the reach of small children!
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect the batteries from excessive heat.
- **⚠ Risk of explosion!** Do not throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a relatively long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries!
- Do not disassemble, open or crush the batteries.

General notes

- Please note that measuring tolerances are possible for technical reasons, as this scale is not calibrated for use in a professional medical context.
- The scale's maximum capacity is 200 kg (440 lb/31 st). The results for weight measurement are displayed in increments of 50 g (0.1 lb).
- Measurements of body fat, body water and muscle content are displayed in increments of 0.1%.
- The calorie requirement is indicated in increments of 1 kcal.
- When supplied to the customer, the scale is set to weigh and measure in "cm" and "kg". You can change the unit settings during initial use with the app.
- Place the scale on an even, hard surface; a hard surface is crucial for achieving accurate measurements.
- Repairs may only be carried out by Beurer Customer Services or authorised retailers. Before submitting a complaint, please check the batteries first and replace them if necessary.

Storage and maintenance

The accuracy of the measurements and service life of the device depend on its careful handling:

IMPORTANT

- The device should be cleaned from time to time. Do not use any abrasive cleaning products and never submerge the device in water.
- Make sure that no liquids come into contact with the scale. Never submerge the scale in water. Never rinse it in running water.
- Do not place any objects on the scale when it is not in use.
- Do not press the button violently or with pointed objects.
- Do not expose the scale to high temperatures or strong electromagnetic fields (e.g. mobile telephones).
- Protect the device from knocks, damp, dust, chemicals, marked temperature fluctuations and nearby sources of heat (ovens, heaters).

General tips


- As far as possible, weigh yourself at the same time every day (ideally in the morning), after having been to the toilet, on an empty stomach and without clothing, in order to obtain results which can be compared.
- Important point regarding the measurement: the calculation of body fat may only be carried out using measurement with the hands and bare feet; the soles of the feet may be lightly dampened where appropriate.
- Results may be unsatisfactory if the soles of the feet are completely dry or have a large amount of hard skin since conductivity will be impaired.
- Stand up straight and still during measurement.
- Wait a few hours after any physical exertion to which your body is unaccustomed.
- Wait approx. 15 minutes after getting up so that the water stored in your body can disperse.
- Remember that only the long-term trend is important. Short-term changes in weight within a few days are normally caused by loss of fluids; however, body water plays an important role in terms of well-being.

Limitations

When calculating body fat and the other values, unusual and implausible results may occur in the case of:


- Children under approx. 10 years of age
- Persons with a fever, persons undergoing dialysis, persons with oedema symptoms and persons suffering from osteoporosis
- Persons taking cardiovascular medication (affecting the heart and vascular system)
- Persons taking vasodilative or vasoconstrictive medication
- Persons with significant anatomical abnormalities in their legs compared with their overall body height (leg length significantly shorter or longer)
- Performance athletes and bodybuilders.

4. Device description

1. Display
2. > Menu: right button (sensor button)
3. Electrodes, right hand
4. ▾ Menu: down button (sensor button)
5. ▲ Menu: up button (sensor button)
6. SET Set button, confirm the entry (sensor button)
7. ⏻ ON/OFF button (sensor button)
8. Electrodes, left hand
9. < Menu: left button (sensor button)
10. Aluminium stand
11. Foot electrodes, right foot
12. Foot electrodes, left foot

Display


1. Date	9. Calorie consumption
2. Time	10. Body fat percentage with interpretation
3. Bluetooth®	11. Body water content with interpretation
4. General overview	12. Muscle percentage with interpretation
5. Segment analysis	13. Bone mass with interpretation
6. Graphic representation	14. Body weight with BMI interpretation
7. User memory	15. Weight unit (kg, lb or st)
8. Average values	


5. Initial use

Positioning the scale


Place the scale on a flat and stable surface. A solid surface is a prerequisite for a precise measurement.


1. Open the cardboard box.


2. Remove the polystyrene cover on the left and right.


3. Lift up the upper and lower part of the aluminium stand at the same time.


4. Insert the upper part of the aluminium stand into the lower part and secure with screws at the top.


5. Position the aluminium stand upright in the cardboard box. Slide the scale to the edge...

6. ...and secure with screws at the bottom.

7. Hold the stand in place with one hand and remove the film with the other.

8. Insert batteries the correct way round in the underside of the scale (+/-). Stand the scale upright on a flat surface.

6. Configuring the diagnostic bathroom scale with an app

To calculate your body fat percentage and additional body values, you must save the personal user data on the scale.

The scale has 10 user memory spaces which allows, for example, you and your family members to save personal settings.

It is also possible to activate users via other mobile end devices on which the "beurer HealthManager" app has been installed or by changing the user profile in the app (see the app settings).

Keep the smartphone close to the scale in order to maintain an active *Bluetooth*[®] connection during initial use.

- **Activate *Bluetooth*[®] in your smartphone's settings.**
- **Download the free "beurer HealthManager" app from the Apple App Store (iOS) or from Google Play (Android[™]).**
- **Start the app and follow the instructions.**
 - The following settings must be set or entered in the "beurer HealthManager" app:

User data	Set values												
Initials	Max. 3-digit (A-Z, 0-9)												
Gender	Male (♂), female (♀)												
Height	100 to 250 cm (3' 03" to 8' 02")												
Date of birth	Year, month, day												
Activity level	The medium-term and long-term perspective is decisive when selecting the activity level:												
	<table border="1"> <thead> <tr> <th></th> <th>Physical activity</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>None.</td> </tr> <tr> <td>2</td> <td>Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).</td> </tr> <tr> <td>3</td> <td>Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.</td> </tr> <tr> <td>4</td> <td>High: physical exertion, at least 4 to 6 times a week for 30 minutes each.</td> </tr> <tr> <td>5</td> <td>Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.</td> </tr> </tbody> </table>		Physical activity	1	None.	2	Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).	3	Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.	4	High: physical exertion, at least 4 to 6 times a week for 30 minutes each.	5	Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.
		Physical activity											
	1	None.											
	2	Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).											
	3	Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.											
4	High: physical exertion, at least 4 to 6 times a week for 30 minutes each.												
5	Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.												
Target weight	kg/lb/st (for progress graph)												

- The app assigns the next free user memory space (person 1 to person 10) to you.
- Carrying out the measurement: For automatic personal recognition, the first measurement must be assigned to your personal user data. The app will prompt you to do this. For how to do this, please see chapter "8. Carrying out the measurement".

7. Alternatively: configuring the diagnostic bathroom scale without an app

We recommend using the app for the initial use and implementing all settings. See chapter "6. Configuring the diagnostic bathroom scale with an app". However, this can also be carried out on the device itself.

Specify device/time/user settings on the control unit. To do so, use the [] and [] buttons to make your selection and confirm with [SET] respectively.

a.) Setting the time

[SET] → "SET TIME" → [SET]

Sequence	Set values
Date	Year–Month–Day
Time	Hours–Minutes

b.) Specifying device settings

[SET] → "SET UNIT" → [SET]

Sequence	Set values
Language	DE EN FR ES IT
Measurement unit	kg–lb–st
Time	24 h–12 h
Guest mode	ON–OFF
Bluetooth®	ON–OFF  When set to "OFF", a Bluetooth® connection is only possible if the scale is switched on. This can be useful for extending the battery life or to prevent undesirable radio radiation.

c.) Setting user data

[SET] → "SET USER" → [SET]

User data	Set values
Initials	Max. 3-digit (A-Z, 0-9)
Gender	Male (♂), female (♀)
Height	100 to 250 cm (3' -03" to 8' -02")
Date of birth	Year, month, day
Activity level	1 to 5 (see section 7)
Target weight	kg/lb/st (for progress graph)
User PIN	Assigned automatically for Bluetooth® use at a later time. This information is only required for app set-up in the future, if existing scale users are to be added during app installations at a later time.

Logging in to the app at a later time

If you have carried out the initial use of the scale without using the app, you can add your existing scale user to the "beurer HealthManager" app. This requires you to enter your 4-digit user "PIN". This process is described in the app.

8. Carrying out the measurement

8.1 Body weight measurement

- Press the  button on the control unit. 0.0 kg appears in the display.
- The scale is then ready to take a weight measurement. To measure your weight, step onto the scale with your weight distributed evenly on both feet. The weight data will not be stored; user allocation will not be possible, as there is no BIA measurement via the hands or feet.


8.2 Carrying out guest measurement


For one-time measurement or if no user data is to be stored. For this, guest mode must be switched on in the device settings.

 → "Guest" → [SET] or switch scale on


User data	Set values
Gender	Male (♂), female (♀)
Height	100 to 250 cm (3' -03" to 8' -02")
Date of birth	Year, month, day
Activity level	1 to 5 (see section 7)

8.3 User measurement, total body analysis

<ul style="list-style-type: none"> – Press the  button on the control unit. 0.0 kg appears in the display.  If guest mode has been activated in the scale settings, step onto the scale when it is switched off. – Stand on the scale barefoot holding the control unit in your hands. Hold the control unit with the underside of your thumbs on the front electrodes and with your fingers pressed firmly against the back electrodes. – Stand still with your weight distributed equally on both feet. Make sure there is no skin contact between your two feet, legs, calves and thighs. 	
<ul style="list-style-type: none"> – If the scale recognises you automatically, the impedance for the analysis will be measured after the weight measurement. This is indicated by the symbol showing current passing through the body. – As soon as the total body analysis is complete, the overview will be shown. The measurement is automatically stored. 	


-  Note: The measurement will not be saved if the following applies:
- The user is not automatically detected. This happens if your measured weight differs from the previous measurement by +/- 3 kg or if another user has similar values (+/- 3 kg) to yours. In this case, press the  button to select the user manually and to carry out the measurement again.
 - The user is standing on the scale with shoes on.
 - The user is not holding the handheld unit with their hands.

9. Analysis of measured values






The results will be analysed automatically after taking the measurement or after selecting the user manually on the control unit. Select your user using the  button. Confirm using the [SET] button or wait 3 seconds.

General overview



<p>The general overview appears at a glance in the display either before or after a total body analysis for your user.</p>	
--	--



Segment analysis for body fat and muscle percentage


<p>a) Body fat in left and right arm</p> <ul style="list-style-type: none"> – Press the menu button  until  is selected. The segment analysis for body fat in the left and right arm appears in the display. 	
--	--


<p>b) Body fat in left and right leg</p> <p>– Use the [] button to select the detailed display of body fat for the left and right leg.</p>	
<p>c) Body fat in the torso</p> <p>– Use the [] button to select the detailed display of body fat in the torso.</p>	
<p>d) Visceral fat (organ fat/abdominal fat)</p> <p>– Use the [] button to select the detailed display of visceral fat.</p>	
<p>a) Muscle percentage in left and right arm</p> <p>– Use the [] button to select the detailed display of muscle percentage in the left and right arm.</p>	
<p>b) Muscle percentage in left and right leg</p> <p>– Use the [] button to select the detailed display of muscle percentage in the left and right leg.</p>	
<p>g) Muscle percentage in torso</p> <p>– Use the [] button to select the detailed display of muscle percentage in the torso.</p>	


Graphic representation

<p>a) Body weight</p> <p>– Press the menu button [] until  is selected. The trend for your body weight appears in the display. The trend for your body weight is shown from measurement to measurement. Later measurements are added on the right. The reference line shows your target weight.</p> <p> You can change your target weight in the user data at any time.</p>	
<p>b) Body fat </p> <p>– Use the [] button to switch to showing body fat values .</p> <p>The trend for your body fat values is shown from measurement to measurement. The reference line represents the average value for all body fat measurements; the bars show whether the values have gone up or down.</p>	

<p>c) Water proportion ☹</p> <p>– Use the [☑] button to switch to showing water proportion ☹.</p> <p>The trend for your body's water proportion is shown from measurement to measurement. The reference line represents the average value for all water proportion measurements; the bars show whether the values have gone up or down.</p>	 <p>The graph shows a trend for water proportion over time. The y-axis ranges from 62.3 to 64.3. A horizontal reference line is at 63.3%. The current value is 63.3%. The date is 23-05-2018 and the time is 18:33.</p>
<p>d) Muscle percentage ↘</p> <p>– Use the [☑] button to switch to showing muscle percentage ↘.</p> <p>The trend for your muscle percentage is shown from measurement to measurement. The reference line represents the average value for all muscle percentage measurements; the bars show whether the values have gone up or down.</p>	 <p>The graph shows a trend for muscle percentage over time. The y-axis ranges from 40.4 to 42.4. A horizontal reference line is at 41.4%. The current value is 41.4%. The date is 23-05-2018 and the time is 18:33.</p>

<p>M Showing stored measurements</p> <p>– Press the menu button [➤] until M is selected. The most recent measured stored measurement appears in the display. The measurement date and the memory space number are shown at the top, e.g. (03-03) means the third memory space = the last of a total of 3 occupied memory spaces.</p> <p>– Use the [↶, ↷] buttons to cycle through the memory spaces.</p>	 <p>The screen shows the most recent stored measurement: 75.4 kg. It also displays three other measurements: 65.2%, 66.8%, and 63.4%. The date is 23-05-2018 and the time is 18:33. The memory space number is 03-03.</p>
--	---

<p>☉ Average values</p> <p>– Press the menu button [➤] until ☉ is selected. All average values in the stored measurements are shown in the display.</p> <p>i The dates of the first and the last measurement are shown at the top.</p>	 <p>The screen shows average values: 25.2%, 26.8%, 63.3%, and 41.6%. The date is 23-05-2018 and the time is 18:33. The current value is 75.4 kg.</p>
---	--

<p>kcal Calorie consumption</p> <p>– Press the menu button [➤] until kcal is selected. The theoretically calculated basal metabolic rate BMR and the active metabolic rate AMR are shown in Kcal.</p>	 <p>The screen shows BMR: 1499 KCAL and AMR: 2248 KCAL. The date is 23-05-2018 and the time is 18:33.</p>
---	---

Switching off the scale (standby mode)

– Wait until the display switches itself off or press the [⏻] button on the control unit to switch the display off manually. The date and time remain visible.

10. Evaluating the results

Body fat percentage

Body fat values are shown in the display and analysed with the bar chart. The subsequent body fat percentages serve as a guideline (for further information ask your doctor).

Male

Age	Low	Normal	High	Very high
10-14	<11 %	11-16 %	16.1-21 %	>21 %
15-19	<12 %	12-17 %	17.1-22 %	>22 %
20-29	<13 %	13-18 %	18.1-23 %	>23 %
30-39	<14 %	14-19 %	19.1-24 %	>24 %
40-49	<15 %	15-20 %	20.1-25 %	>25 %
50-59	<16 %	16-21 %	21.1-26 %	>26 %
60-69	<17 %	17-22 %	22.1-27 %	>27 %
70-100	<18 %	18-23 %	23.1-28 %	>28 %

Female

Age	Low	Normal	High	Very high
10-14	<16 %	16-21 %	21.1-26 %	>26 %
15-19	<17 %	17-22 %	22.1-27 %	>27 %
20-29	<18 %	18-23 %	23.1-28 %	>28 %
30-39	<19 %	19-24 %	24.1-29 %	>29 %
40-49	<20 %	20-25 %	25.1-30 %	>30 %
50-59	<21 %	21-26 %	26.1-31 %	>31 %
60-69	<22 %	22-27 %	27.1-32 %	>32 %
70-100	<23 %	23-28 %	28.1-33 %	>33 %

Values will often be lower for fitness enthusiasts. Depending on the sporting activities performed, the training intensity and the physical constitution, the results may still be below the stated standard values. Please note, however, that there may be a risk to health in the case of extremely low values.

Visceral fat (organ fat/abdominal fat)

Visceral fat denotes fat that surrounds the internal organs within the abdominal cavity. Visceral fat is suspected of contributing to conditions such as high blood pressure or diabetes.

Visceral fat is shown in the display as a value between 1 and 59.

Visceral fat percentage value	Rating
1-12	Your visceral fat percentage is in the normal range.
13-59	Your visceral fat percentage is too high.

Body water content

The body water content in % is normally within the following ranges:

Male

Age	Poor	Good	Very good
10-100	<50 %	50-65 %	>65 %

Female

Age	Poor	Good	Very good
10-100	<45 %	45-60 %	>60 %

Body fat contains relatively little water, therefore body water content may be below the standard values in persons with a high body fat percentage. On the other hand, endurance athletes may exceed the standard values due to having low levels of fat and a high muscle percentage.

The body water calculation performed using this scale is not suitable for drawing medical conclusions, for example in the case of age-related water retention. Consult your doctor where necessary. A high body water content is generally desirable.

Muscle percentage

The muscle percentage measurement measures the skeletal muscle mass. No organ muscles are taken into account. The muscle percentage is normally within the following ranges:

Male

			
Age	Low	Normal	High
10-14	<44%	44–57 %	>57 %
15-19	<43%	43–56 %	>56 %
20-29	<42%	42–54 %	>54 %
30-39	<41%	41–52 %	>52 %
40-49	<40%	40–50 %	>50 %
50-59	<39%	39–48 %	>48 %
60-69	<38%	38–47 %	>47 %
70-100	<37%	37–46 %	>46 %

Female

			
Age	Low	Normal	High
10-14	<36%	36–43 %	>43 %
15-19	<35%	35–41 %	>41 %
20-29	<34%	34–39 %	>39 %
30-39	<33%	33–38 %	>38 %
40-49	<31%	31–36 %	>36 %
50-59	<29%	29–34 %	>34 %
60-69	<28%	28–33 %	>33 %
70-100	<27%	27–32 %	>32 %

Bone mass

Like the rest of our body, our bones are subject to natural growth, shrinking and ageing processes. Bone mass increases rapidly during childhood and reaches its peak between the ages of 30 and 40. As we age, our bone mass then begins to reduce. You can combat this reduction to an extent with the help of a healthy diet (particularly calcium and vitamin D) and regular physical exercise. The stability of the skeleton can be further increased by means of targeted muscle build-up. Note that this scale identifies the mineral content of the bone. It is very difficult to influence the bone mass, although it does fluctuate within the scope of influencing factors (weight, height, age, gender). There are no recognised guidelines or recommendations.

Important

Please do not confuse bone mass with bone density. The bone density can only be determined by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions on changes to bones and bone hardness (e.g. osteoporosis) using this scale.

BMR

The basal metabolic rate (BMR) is the amount of energy required by the body at complete rest in order to maintain its basic functions (e.g. lying in bed 24 hours a day). This value is primarily dependent on weight, height and age. It is displayed on the diagnostic scale in the unit kcal/day and calculated using the scientifically recognised Mifflin St. Jeor equation.

This is the amount of energy that is required by your body under all circumstances and must be re-supplied to the body in the form of food. If your energy intake is below this level in the longer term, this can result in damage to health.

AMR

The active metabolic rate (AMR) is the amount of energy that the active body consumes each day. An individual's energy requirement increases when the level of physical activity increases; the diagnostic bathroom scale calculates this by means of the activity level (1–5) that has been entered.

To maintain the current weight, the energy that the body uses must be replaced in the form of food and drink. If less energy is taken in than is used over a long period of time, the body primarily takes the difference from the fat stores with resulting weight loss. However, if the amount of energy taken in exceeds the calculated active metabolic rate (AMR) for a longer period, the body cannot burn off the excess energy. The excess is stored in the body as fat, leading to weight gain.

Temporal context of the results

 Please note that only the long-term trend is significant. Brief deviations in weight within a few days are normally caused by loss of fluids.

The interpretation of the results is based on changes in total body weight, percentage of body fat, body water and muscle content, as well as on the length of time over which these changes occur. Rapid changes within the scope of a few days are to be considered separately from medium-term changes (in the scope of weeks) and long-term changes (months).

It can be said as a basic rule that short-term changes in weight are almost entirely changes in water content, whereas medium-term and long-term changes may also relate to the fat percentage and the muscle percentage.

- If your weight drops in the short term but your body fat percentage increases or stays the same, you have only lost water, for example following a training session, visit to the sauna or a diet aimed only at fast weight loss.
- On the other hand, if your weight increases in the medium term and your body fat percentage drops or stays the same, you may have built up valuable muscle mass.
- If your weight and body fat percentage both fall at the same time, then your diet is working – you are losing fat mass.
- Ideally, you should support your diet with physical activity, fitness or strength training. This enables you to increase your muscle percentage in the medium term.
- Body fat, body water and muscle percentages should not be totalled (muscle tissue also contains components made of body water).


11. Transferring measurements via app

If the app is not open, the newly assigned measurements are saved on the scale. A maximum of 30 measurements per user can be saved on the scale. The saved measurements are transferred automatically to the app when you open the app within the *Bluetooth*[®] range. It is possible to automatically transfer data when the scale is switched off if *Bluetooth*[®] is enabled in the device settings.

12. Additional information


Deleting scale data

a) Delete user with all settings and saved values


- After switching on, select your user by pressing the [] button. If necessary, press the button several times and confirm with [SET].
- When  appears in the display, hold down the [SET] button for 12 seconds to delete the user.
- Use the [, ,] buttons to select “YES” and confirm the prompt using the [SET] button. If you do not want to delete the user, confirm “NO” using the [SET] button.

 All saved values and settings for this user will be deleted.

b) Delete individual saved value for a user

- Select a measured value in the user memory  and briefly press the [SET] button on the control unit to delete the stored value displayed.
- You can also do this immediately after the measurement is taken.
- Use the [, ,] buttons to select “YES” and confirm the prompt using the [SET] button. If you do not want to delete the value, confirm “NO” using the [SET] button.

c) Delete all stored values for a user

- Open the user memory  and press the [SET] button on the control unit for around 5 seconds to delete all values.
- Use the [, ,] buttons to select “YES” and confirm the prompt using the [SET] button. If you do not want to delete, confirm “NO” using the [SET] button.

Changing the batteries

The control unit has a low battery indicator. “LOW BATTERY” appears in the display if the batteries in the scale are too weak. In this case, the batteries should be replaced. Connect the scale to the app again. The system time is updated.

**NOTE:**

- When replacing the batteries, use batteries of the same type, make and capacity.
- Do not use rechargeable batteries.
- Use batteries free from heavy metals.

13. Cleaning and maintaining the device

The device should be cleaned from time to time.

This should be done using a damp cloth and, if necessary, a small amount of detergent.

**IMPORTANT**

- Never use harsh solvents or cleaning products!
- Never submerge the device in water!
- Do not clean the device in the dishwasher.

14. Disposal

Empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

For environmental reasons, do not dispose of the device in the household waste at the end of its service life. Dispose of the device at a suitable local collection or recycling point in your country.

Dispose of the device in accordance with EC Directive–WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Dispose of packaging in an environmentally friendly manner.

**15. What if there are problems?**

If the scale encounters an error during measurement, the following is displayed:

Display	Cause	Solution
ERROR:_OVERLOAD	Maximum weight capacity of 200 kg exceeded.	Loads must be limited to 200 kg.
LOW_BATTERY	Batteries are nearly flat.	Replace the batteries (see section on changing the batteries).
MEMORY_FULL	30 measurements have been stored.	Transfer the measured values.
ERROR:_CONTACT ERROR	The electrical resistance between the electrodes and the soles of the feet is too high (e.g. with very hard skin).	Repeat the measurement barefoot. Dampen the soles of your feet slightly. If necessary, remove hard skin on the soles. Hold the hand electrodes with the fingers behind and the balls of the thumb in front.
ERROR:_DATA	The BIA measurement was successful, but the fat percentage is outside the measurable range (less than 3% or greater than 65%).	Please repeat the measurement barefoot, or with the soles of the feet lightly dampened if necessary. If necessary, change your set activity level.
Empty display	Batteries are completely flat.	Replace the batteries (see section on changing the batteries).
	The batteries are not inserted properly.	Check for correct polarity, (see "Initial use" -> "Changing the batteries").

Display	Cause	Solution
Only weight is displayed, no diagnosis	Unknown measurement as it is outside the user assignment limit or a unique assignment is not possible.	Select user on the scale using the [] button or repeat user assignment in the app. Hold the electrodes of the handheld unit and stand barefoot on the scales.

Incorrect weight measurement

Possible faults	Solution
Incorrect zero point of scale.	Switch the scale on, wait approx. 5 seconds and repeat the measurement.
Scale is on a carpet.	Place the scale on a flat, solid surface. Stand on the outside of the scale as far as possible without it tipping.
Incorrect weight is displayed.	Stand as still as possible.
Incorrect weight is displayed.	Wait until the scale switches itself off. Then repeat the measurement.

No *Bluetooth*[®] connection (ⓧ symbol is missing)

Possible faults	Solution
Radio connection	The range in the open is approx. 10 m. Walls and ceilings reduce the range. Other radio waves may interfere with the transmission. Therefore, do not position the scale near devices such as WIFI routers, microwaves, induction hobs etc.
<i>Bluetooth</i> [®] connection	Completely close the app (including in the background). Switch <i>Bluetooth</i> [®] off and back on again. Switch the smartphone off and back on again. Remove the batteries in the scale for a short period and reinsert them. See FAQs at www.beurer.com .
Connection set-up taking too long	Faster <i>Bluetooth</i> [®] interval when the scale is switched on.

16. Technical specifications

Scale dimensions	380 x 390 x 1000 mm
Scale batteries	4 x 1.5V type AA
Measurement range	8 kg to 200 kg
Scale interval d:	The display can be read in increments of 0.05 kg.
Repetition accuracy:	The measuring tolerance for repeated measurements is +/- 0.3 kg (several measurements in direct succession on the same scale with position of scale and person the same as far as possible).
Absolute precision:	In comparison with a calibrated weight, the measured value is +/-0.5% +0.05 kg. E.g. at 50 kg this corresponds to +/- 0.3 kg; at 100 kg this corresponds to +/-0.55 kg.
Data transfer via <i>Bluetooth</i> [®] wireless technology	This product uses <i>Bluetooth</i> [®] low energy technology, frequency band 2.402–2.480 GHz, transmission power max. 0 dBm, compatible with <i>Bluetooth</i> [®] 4.0 smartphones/tablets

A list of compatible smartphones, plus information and software for the devices can be found at the following link:
connect.beurer.com

Subject to technical changes.

We hereby confirm that this product complies with the European RED Directive 2014/53/EU.

The CE Declaration of Conformity for this product can be found under:

<https://www.beurer.com/web/we-landingpages/de/cedeclarationofconformity.php>

17. Warranty / service

Beurer GmbH, Söflinger Straße 218, 89077 Ulm, Germany (hereinafter referred to as "Beurer") provides a warranty for this product, subject to the requirements below and to the extent described as follows.

The warranty conditions below shall not affect the seller's statutory warranty obligations which ensue from the sales agreement with the buyer.

The warranty shall apply without prejudice to any mandatory statutory provisions on liability.

Beurer guarantees the perfect functionality and completeness of this product.

The worldwide warranty period is 5 years, commencing from the purchase of the new, unused product from the seller.

The warranty only applies to products purchased by the buyer as a consumer and used exclusively for personal purposes in the context of domestic use.

German law shall apply.

During the warranty period, should this product prove to be incomplete or defective in functionality in accordance with the following provisions, Beurer shall carry out a repair or a replacement delivery free of charge, in accordance with these warranty conditions.

If the buyer wishes to make a warranty claim, they should approach their local retailer in the first instance: see the attached "International Service" list of service addresses.

The buyer will then receive further information about the processing of the warranty claim, e.g. where they can send the product and what documentation is required.

A warranty claim shall only be considered if the buyer can provide Beurer, or an authorised Beurer partner, with

- a copy of the invoice/purchase receipt, and
- the original product.

The following are explicitly excluded from this warranty:

- deterioration due to normal use or consumption of the product;
- accessories supplied with this product which are worn out or used up through proper use (e.g. batteries, rechargeable batteries, cuffs, seals, electrodes, light sources, attachments and nebuliser accessories);
- products that are used, cleaned, stored or maintained improperly and/or contrary to the provisions of the instructions for use, as well as products that have been opened, repaired or modified by the buyer or by a service centre not authorised by Beurer;
- damage that arises during transport between manufacturer and customer, or between service centre and customer;
- products purchased as seconds or as used goods;
- consequential damage arising from a fault in this product (however, in this case, claims may exist arising from product liability or other compulsory statutory liability provisions).

Repairs or an exchange in full do not extend the warranty period under any circumstances.

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The *Bluetooth*[®] word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.