Therabody

RecoveryThermTM

Hot Vibration Back & Core

User Manual

Table of contents

EN	2 - 6	ES	12 - 16	DE	22 - 26
FR	7 - 11	IT	17 - 21	SC	27 - 31

Product Overview

RecoveryTherm™

Hot Vibration Back & Core

Heat and vibration for advanced relief.

From the creators of Theragun® comes RecoveryTherm™ Hot Vibration Back & Core, precisely delivering heat and vibration therapy in one wearable garment. The RecoveryTherm™ Hot Vibration Back & Core is an innovative, battery powered heating technology that uses carbon-fiber mesh and germanium fabric to provide 360 degrees of targeted, consistent heat.

Intended Use


The RecoveryTherm™ Hot Vibration Back & Core is intended to offer a pain management solution for everyone, from alleviating lower back pain to reducing tension or providing relief from menstrual pain. The precise and targeted heat in tandem with vibration therapy helps relieve lower back pain and muscle stiffness and improves blood flow all in one garment.

Getting started

Step 1

Put on your garment

1. Find a comfortable position sitting down or standing upright; we recommend standing for optimal placement.
2. To put the RecoveryTherm™ Hot Vibration Back & Core on, center the Therabody logo across the treatment area, e.g., lower back.
3. Hold each side of the straps with the control module in your right hand. Adjust gently to ensure they are approximately the same distance and that the device is centered.
4. Once centered, wrap the left side of the strap around your abdomen toward your hip and secure it by placing the right strap on top. Secure in place with the Velcro strap. The garment should feel snug but not too tight. The control module should rest between the center of your body and your right hip.

Getting started

Step 2

Use your garment

1. Insert the battery into the control panel. Ensure that the Therabody logo faces outwards and the charging port is on the bottom. Battery insertion should be seamless and not forced.
2. Press and hold the power button for two seconds to turn on the device.
3. RecoveryTherm™ Hot Vibration Back & Core delivers heat on the front and back of your core for a customizable 360-degree heat therapy treatment. Short press the heat buttons to toggle between the three heat settings on both the front and back. The bottom heat button controls the front, and the top heat button controls the back.
4. Add vibration to your back treatment by pressing the vibration button. Toggle between the three vibration modes (low intensity, high intensity, and wave) by short-pressing the vibration button.
5. The treatment will automatically shut off after approximately 20 minutes. Turn the garment back on if you'd like to repeat the treatment.

Step 3

Feel the difference.® Every day.

Whether used for targeted pain management, muscle recovery, or general soreness, RecoveryTherm™ Hot Vibration Back & Core will help keep you feeling your best.


See Warnings for important safety instructions

Scan the QR code to get started.

Getting to know your wearable device

- A** Wrap garment
- B** Removable battery
- C** USB-C Port
- D** Travel neoprene pouch

Modes Overview

Mode 1 - Front heat: There are three temperatures for front heat: low, medium, and high. Short press to toggle between the three temperatures.

Mode 2 - Back heat: There are three temperatures for back heat: low, medium, and high. Short press to toggle between the three temperatures.

Mode 3 - Vibration: There are three intensities for vibration: low, high, and wave. Short press to toggle between the three intensities.


Charging

- To charge the RecoveryTherm™ Hot Vibration Back & Core, remove the battery from the control panel and charge using the provided cable.
- The RecoveryTherm™ Hot Vibration Back & Core includes a Therabody USB-C to USB-A cable.
- The RecoveryTherm™ Hot Vibration Back & Core is USB-C enabled.
- The battery level is displayed in the power button and on the LED indicator on the battery pack, next to the USB-C connector:
 - **Orange: low charge, blue: mid charge, green: fully charged.**
 - **When charging, the LED will flash to indicate it is properly connected and will turn solid green when complete.**

Note: Ensure the charger is from a certified manufacturer and has not suffered any structural damage. Do not try to use the garment while charging.

Cleaning Instructions:

1. The RecoveryTherm™ Hot Vibration Back & Core is NOT waterproof or machine washable. Do not allow excess liquids or water to enter the product.
2. Remove the battery before cleaning.
3. To clean, wipe the inner surface of the strap with a dampened microfiber cloth.

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